

With Newlight adaptor ADA-302, a waveplate with OD=30.0mm can be mounted to EITHER  
(a) Thorlabs RSP1 - Rotation Mount for  $\varnothing$ 1" Optics, with the SM1 thread, OR  
(b) Thorlabs RSP2 - Rotation Mount for  $\varnothing$ 2" Optics, with a 2" retaining ring

for rotation of the principal axis of the waveplate.



Adaptor ADA-302 with an OD 30.0mm waveplate mounted

### Thorlabs RSP1




Thorlabs SM1 (1.035"-40) thread



ADA-302

Dia=50mm



### Thorlabs RSP2